

FEATURES

NO Heavy Tanks NO heavy Carts

- Continuous oxygen supply
- Pulse flow operation
- · Battery operated
- Purity: 93%+ -3% at all settings
- Light weight: 1.50kg single battery
- Light weight: 1.76kg double Battery
- Light weight: 1.99kg triple Battery
- Dimensions:
 - Length: 157mmWidth: 86mmHeight: 170mm
- Smart Intelligence If the POC cannot detect breathing, it will automatically start to pulse every 3 seconds.
- Synchronised ability- adjust according to patients breathing pattern.
- 5l Flow settings





PSA TECHNOLOGY



O2 CAPACITY TO **1**L



1.5KG LIGHT WEIGHT



7.5 HOURS
BATTERY RUN TIME with a tripple battery

1.50 KGSmallest and Lightest1 L-Oxygen concentrator in the world.

Equivalent to 5 Lites, Designed for 24/7 use.

5cm large LCD color display screen easy to read.




FLOW SETTINGS



Flow settings 1-5 can meet different oxygen demands.

Setting 1 210 ml/min
Setting 2 420 ml/min
Setting 3 600 ml/min
Setting 4 800 ml/min
Setting 5 1000 ml/min

PRODUCT INTRODUCTION



EXTERNAL BATTERY CHARGER

Model # BC-K

when the battery is set correctly , the flashing blue LED will indicate that battery is charging.



when the green light illuminated, the battery is fully charged.



BATTERY DURATION

*Battery	Single Battery	Double Battery	Tripple Battery
Duration	SKU# BA-K201	SKU# BA-K200	SKU# BA-K202
Flow Setting Battery duration	1: 2.75 h 2: 2.00 h 3: 1.50 h 4: 1.30 h 5: 1.10 h	1: 5.50 h 2: 3.75 h 3: 3.00 h 4: 2.50 h 5: 2.20 h	1: 7.05 h 2: 5.20 h 3: 4.00 h 4: 3.50 h 5: 2.80 h

^{*}Based on 20 BPM. Actual duration results may vary.

Battery duration times are based on new, fully charged batteries. Duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the portable concentrator.

ACCESSORIES

Standard package:

Carry Bag, Nasal Cannula, AC & DC Adapter, Single Battery, Filters



Optional Accessories:

External battery charger, Sieve beds, Double Battery, Triple Battery

